

STOCCARDA








mattress with **ZONES** of different bearing capacities





- 3,5 cm side thermosensitive memory form, self-modelling to distribute the body weight evenly.
- **B** Perimeter casing against sinking, perforated to allow a constant and healthy air recirculation.
- **G** Spring system with 1600 independent pocketed micro-springs with 7 zones bearing differentiated capacity.



• High density embossed layers with over 1000 pressure points for optimal comfort.


















Mattress Maintenance:

Periodically rotate head / feet, up / down and ventilate.









Elaborate side band complete with 3D breathable insert and fabric handles for the mattress positioning on the bed frame.



STOCCARDA	190 / 195 / 200						
	80	85-90	120	140	160	165-170-180	per cm²
Silver-Tencel - Cashmere - Bamboo - Hypoallergenic - Organic Cotton							
Easy							

Warning: the medical device is NOT produced with the Easy lining.

Mattress Maintenance:

Periodically rotate head / feet, up / down and ventilate.



Fabric with silver ion treatment for antibacterial and odour protection. Under-edge perimeter zip for complete removing of linings.



Elaborate side band complete with 3D breathable insert and fabric handles for the mattress positioning on the bed frame.



STOCCARDA ELITE		Taylor made					
	80	85-90	120	140	160	165-170-180	per cm²
Silver-Tencel – Cashmere – Bamboo – Hypoallergenic – Organic Cotton							
Base							

